



INTRAORAL SCANNER SETUP INSTRUCTIONS

iTero Element

1. Call iTero support (1-800-577-8767) and ask them to have the scans sent to OrthoSelect. Give them the company ID, 14950, and they will instruct you further.
2. When starting a scan, select "New Scan", and enter patient information.
 - a. Choose Case Type: "iRecord"
 - b. Choose Ship To: "OrthoSelect"



The screenshot shows the 'New Scan' window in the iTero software. The window title is 'New Scan' and the user is 'Dr. Test, Test'. The license number is '12345'. The form contains the following fields: 'First Name', 'Last Name', 'Chart #', 'Case Type' (a dropdown menu), 'Due Date', and 'Ship To' (a dropdown menu). A red box highlights the 'Case Type' and 'Ship To' dropdown menus. Below the form is a 'Notes' section with an 'Add Notes' button.



3Shape TRIOS

1. Log in to your account at us.3shapecommunicate.com.
2. Click on "Connections".
3. Click on "Add Connection".
4. Type the email address intraoral@myorthoselect.com and click "Find".

CareStream 3600

1. Tell your CareStream representative that you want all of your scans to go through CareStream's system.
2. Provide them with the email address intraoral@myorthoselect.com, which will send a connection invitation to OrthoSelect.

Medit i500

1. Click [here](#) for instructions.

3M True Definition

1. Call 3M support at 1-800-634-2249, select option 3, then option 1.
2. Request to send intraoral scans to OrthoSelect.
3. Provide 3M representative with our location, American Fork, UT, and email address, intraoral@myorthoselect.com.
4. When starting a scan, select a date request, then select “OrthoSelect”.

Other Scanners

1. For most other scanners, you should be able to get separate upper and lower arch files, which need to be in .stl format.
2. If possible, make sure files are oriented in occlusion.
3. Use “Direct Upload” option on the Rx form to upload .stl files.