

Recommended Products

ETCHANT



Blue Gel Etchant by Reliance

- Controlled Application
- 35% Phosphoric Acid
- 800-323-4348
- relianceorthodontics.com

PRIMER



Assure Plus by Reliance

- Strong Bonds on All Surfaces
- Compatible with Etch Prep
- 800-323-4348
- relianceorthodontics.com

ADHESIVE



L.C.R. or Flow Tain by Reliance

- Flowable Consistency with the Strength of a Composite
- Hybrid Flowable
- 800-323-4348
- relianceorthodontics.com

UV CURING LIGHT



Ledex WL-090 by Practicon

- Precise Power and Duration Cures the Adhesive Fast.
- 800-959-9505
- practicon.com

Alternatively, you may scan the following QR Code and enter in the item numbers to add them to your cart. These are the products we recommend using when bonding with DIBS AI trays.

5g Blue Gel Etchant Syringe w/ Tips.....#BEGS

6ml Assure Plus All Surface Bonding Resin.....#PLUS

0.2ml Assure Plus 50ct Unidose (10ml).....#UAP50

5g L.C.R. Paste in Push Syringe.....#LCR

1.5g L.C.R. in Push Syringe w/ Tips.....#LCR15

5g Flow Tain Paste in Push Syringe.....#FT

